

Tel: +44 (0) 1306 884885 Fax: +44 (0) 1306 884886

FrameStar® 480/384

Description: 384 well PCR plate for the Roche LightCycler® 480. Rigid two-component design eliminates warping and distortion during PCR making it ideal for use with robotic systems.

Ordering Information:	4ti-0380/C	FrameStar® 480/384 clear frame, clear wells
	4ti-0381	FrameStar® 480/384 clear frame, white wells
	4ti-0382	4ti-0381 & 50 qPCR adhesive seals (4ti-0560)
	4ti-0383	4ti-0380/C & 50 qPCR adhesive seals (4ti-0560)

Sealing Methods FrameStar® 480/384 can be sealed using the following seals:

Adhesive

4ti-0500	PCR Seal, 100 sheets (135mm x 80mm)
4ti-0501	PCR Seal, 1 roll (100m x 80mm)
4ti-0512	Peelable Solvent Resistant Foil, 100 sheets (122mm x 80mm)
4ti-0515	Gas Permeable Seal, 50 sheets (135mm x 80mm)
4ti-0515/ST	Gas Permeable Seal, 50 sheets (135mm x 80mm), sterile
4ti-0515/R	Gas Permeable Seal, 1 roll (100m x 80mm)
4ti-0550	PCR Foil Seal, 100 sheets (135mm x 80mm)
4ti-0551	PCR Foil Seal, 1 roll (200m x 80mm)
4ti-0560*	qPCR Seal, 100 sheets (135mm x 80mm)
4ti-0561*	qPCR Seal, 1 roll (100m x 80mm)
4ti-0565*	Q-Stick™ for qPCR, 100 sheets, (135mm x 80mm)

Heat

4ti-0520	Peel Seal 600m x 78mm roll
4ti-0521	Peel Seal 500m x 115mm roll
4ti-0522	Peel Seal 100 sheets (125mm x 78mm)
180°C for 2 seconds	
4ti-0530	Pierce Seal 610m x 78mm roll
4ti-0530/REMP	Pierce Seal 610 x 78mm roll
4ti-0531	Pierce Seal 100 sheets (125mm x 78mm)
4ti-0532	Pierce Seal 500m x 115mm roll
4ti-0532/REMP	Pierce Seal 500m x 115mm roll
175°C for 2 seconds	

Tel: +44 (0) 1306 884885 Fax: +44 (0) 1306 884886

4ti-0535	Foil Seal 610m x 78mm roll
4ti-0535/REMP	Foil Seal 610m x 78mm roll
4ti-0536	Foil Seal 100 sheets (125mm x 78mm)
4ti-0537	Foil Seal 500m x 115mm roll
4ti-0537/REMP	Foil Seal 610m x 78mm roll
175°C for 2 seconds	
4ti-0540*	Clear Seal 500m x 78mm roll
4ti-0540/REMP*	Clear Seal 500m x 78mm roll
4ti-0542*	Clear Seal 350m x 115mm roll
4ti-0542/REMP*	Clear Seal 500m x 115mm roll
4ti-0541*	Clear Seal 100 sheets (125mm x 78mm)
180°C for 2 seconds	
4ti-0573*	Clear Weld Seal Mark II 610m x 78mm roll
4ti-0574*	Clear Weld Seal Mark II 500m x 115mm roll
4ti-0575*	Clear Weld Seal Mark II 100 sheets (125mm x 78mm)
180°C for 2 seconds	
4ti-0580	Clear Seal 3730 610m x 78mm roll
4ti-0582	Clear Seal 3730 500m x 115mm roll
4ti-0581	Clear Seal 3730 roll of 1000 perforated sheets
180°C for 2 seconds	
4ti-0585	DMSO Resistant Peel Seal 500m x 78mm roll
4ti-0585/REMP	DMSO Resistant Peel Seal 500m x 78mm roll
4ti-0586	DMSO Resistant Peel Seal 500m x 115mm roll
4ti-0586/REMP	DMSO Resistant Peel Seal 500m x 115mm roll
4ti-0587	DMSO Resistant Peel Seal 100 sheets (125mm x 78mm)
180°C for 2 seconds	
4ti-0590	Thermal Bond 500m x 78mm roll
4ti-0591	Thermal Bond 100 sheets (125mm x 78mm)
4ti-0592	Thermal Bond 200m x 115mm roll
180°C for 2 seconds	
4ti-0596	Gas Permeable Seal 200m x 78mm roll
4ti-0595	Gas Permeable Seal 100 sheets (125mm x 78mm)
175°C for 2 seconds	

Storage Conditions: Room temperature. Avoid direct exposure to light and high humidity.

Material Safety Data: Non-hazardous