



The new Tear-A-Way™ plate is perforated for easy division into strips. The non-skirted universal PCR plate of choice.

Tear-A-Way™ plates allow for the most flexible and efficient use of a PCR plate. No need to run half-empty plates or to cut plates with scissors. Cutting plates risks perforating wells, damaging sealing rings and contamination.

Tear-A-Way™ plates can be quickly and easily separated along the perforations between the rows. Every single strip can be utilised, saving time and costs. Please see our website to view the compatibility of Tear-A-Way™ with your thermal cycler or sequencer.



Features and Benefits

- Plate is perforated to enable accurate tearing into strips
- Numbers on all strips to prevent sample errors
- Non-skirted universal cycler and sequencer compatibility
- Sealable with 4titude PCR cap strips (4ti-0751 and 4ti-0752), adhesive seals and heat seals
- Simply snaps to fit a 24 or 48-well thermal cycler block

Code	Description	Quantity
4ti-0750/TA	Tear-A-Way™ PCR plate, non-skirted, clear	50 plates
4ti-0750/W/TA	Tear-A-Way™ PCR plate, non-skirted, white	50 plates