



KIMBLE® KONTES® Dounce Tissue Grinder, with two glass pestles and tube

- Manufactured from borosilicate glass 3.3
- Designed primarily for cellular work where the nucleus remains intact after homogenization.
- All-glass construction

- Two pestles are supplied with each complete unit
- Large clearance pestle is used for the initial sample reduction
- Small clearance pestle is used to form the final homogenate
- Replacement components are available and completely interchangeable

Recommended for 2 mL Tissue Grinder

DOMINIQUE DUTSCHER SAS