

## Datenblatt - *Data Sheet*

*Article* Tube with thread, Round, 20x125 mm, no Cap  
Boro 5.1

*Article Number* 90P20125RW (90405015)

*Material:* **Borosilicate Glass 5.1**

- Water resistance tested to ISO 719 (at 98 °C): HGB 1
- Alkali resistance tested to ISO 695: class A2
- Acid resistance tested to ISO 12116: S1
- Maximum using temperature long-time load 400°C (\*)
- Maximum using temperature Short-time load 500°C (\*)
- Autoclavable

*Dimension/Standard:*

| Over all length H | Outside diameter D | Thread T   | Wall W |
|-------------------|--------------------|------------|--------|
| 125 mm            | 20 mm              | GPI 18-415 | 1,2 mm |

