



SECTION A-A

(\varnothing .49±.01)
(100 PLCS.)

UNLESS OTHERWISE SPECIFIED
DIMENSIONS ARE IN INCHES
TOLERANCES ARE:
FRACTIONS DECIMALS ANGLES
±1/64 .XX ± .01 ± 0.1
.XXX ± .001
.XXXX ± .0005

MATERIAL	APPROVALS	DATE
FINISH	DRAWN	9/28/2005
DO NOT SCALE DRAWING	CHECKED	
	RESP ENG	
	MFG ENG	
	QUAL ENG	

CAD GENERATED DRAWING.
DO NOT MANUALLY UPDATE

APPROVALS	DATE
DRAWN S. KRAMER	9/28/2005
CHECKED	
RESP ENG	
MFG ENG	
QUAL ENG	



1310 Thurman Street
Lodi CA, 95240
Ph: (209) 333-2120
Fx: (209) 333-8623
info@scientificspecialties.com

5540-XX
RACK, 100 PLACE W/HINGED LID

SIZE	SHRINK RATE	REV.
A		-
SCALE 1:2	JOB NO.	SHEET 1 OF 1

PRINTED 12/27/2006