


GCA 358647

THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION AND IS A COPYRIGHT OF WHEATON INDUSTRIES INC. IT MAY NOT BE REPRODUCED OR TRANSMITTED IN ANY FORM OR BY ANY MEANS, ELECTRONIC OR MECHANICAL, INCLUDING PHOTOCOPYING, RECORDING, OR BY ANY INFORMATION STORAGE AND RETRIEVAL SYSTEM, WITHOUT THE EXPRESS WRITTEN CONSENT OF WHEATON INDUSTRIES INC. ALL RIGHTS RESERVED.



THIRD ANGLE PROJECTION



PRODUCT COMPONENTS

ITEM	QTY	DESCRIPTION
1.	144	W012027 TUBE E/B
2.	144	W000036 CAP B/A

UNCONTROLLED COPY



- △ TUBE MATERIAL: TYPE I CLEAR BOROSILICATE GLASS
- CAP MATERIAL: BLACK PHENOLIC WITH .005
- ◇ PTFE FACED 14B WHITE RUBBER LINER

G	2011101211223320	10/31/11 TP
F	2010040911324290	04/20/10 BL
E	2010011410081214	01-21-10 TRJ
D	14145	11-08-02 BL
C	12078	12-05-00 DB
B	10599	08-27-99 DB
A	4822	10-12-94 TRJ
REV	E.C.N. #	DATE

DATE DRAWN 10-26-89 MAT'L  
 DRAWN BY THOMAS R. JOHNSON SL APPROVE. WRT.  
 CHECKED BY THOM LANE AUTHORIZED R.W.S.

TITLE 16 X 125mm TUBE, CULTURE, 15-415, PTFE  
 CAPACITY SIZE 16 X 125mm  
 FINISH G.P.I. 15-415

WHEATON Industries Inc.  
 1501 North 10th Street, Verona, NJ 08882-2038

DWG. NO: B- 358647

SCALE: FULL UNLESS NOTED. DO NOT SCALE THIS DRAWING.

-----NOTES-----