



:\GCA\W216832

COPYRIGHT BY 2008
WHEATON INDUSTRIES INC. ALL RIGHTS
RESERVED. THIS MATERIAL MAY NOT
BE REPRODUCED FOR ANY PURPOSE
WITHOUT THE EXPRESS WRITTEN
CONSENT OF WHEATON INDUSTRIES INC.

THIRD ANGLE PROJECTION

PRODUCT COMPONENTS		
ITEM	QTY	DESCRIPTION



UNCONTROLLED COPY

MATERIAL: TYPE III CLEAR SODA-LIME GLASS

- 2. DIMENSIONS IN [BRACKETS] ARE MILLIMETERS.
- 1. 123cc APPROX BON.

TITLE 4oz BTL, BS RD, TYPE 3 CLR, W/O CAP, BULK	DATE DRAWN 10/17/08 DRAWN BY B. LLOYD	DATE APPROX. VOL. 3.500oz (99.2gr)	REV A E.C.N. #	DATE 10/17/08 DL
CAPACITY 4oz (124.5cc-130.5cc OVERFLOW)	CHECKED BY JB 10/17/08	APPROVED TRJ 10/17/08	REV	DATE
FORM G.P.I. 22-400-405 C/R	SCALE: FULL UNLESS NOTED DO NOT SCALE THIS DRAWING.	TOLERANCE:	DWG. NO: B- W216832	

