



THIS DRAWING IS THE PROPERTY OF WHEATON CORPORATION. IT IS TO BE USED ONLY FOR THE PURPOSES SPECIFIED IN THE ORDER. IT IS NOT TO BE REPRODUCED OR TRANSMITTED IN ANY FORM OR BY ANY MEANS, ELECTRONIC OR MECHANICAL, INCLUDING PHOTOCOPYING, RECORDING, OR BY ANY INFORMATION STORAGE AND RETRIEVAL SYSTEM. WHEATON CORPORATION, 2001052114301198, 06-22-11 TP



THIRD ANGLE PROJECTION



MATERIAL: TYPE III CLEAR SODA-LIME GLASS

2. DEPRESSED THREAD AT SEAM.

1. DIMENSIONS ARE IN INCHES AND [MILLIMETERS].

NOTES

PRODUCT COMPONENTS	
ITEM QTY	DESCRIPTION

TITLE	DATE DRAWN	REVISED	BY	DATE
16 OZ BTL, FRENCH SQUARE, TYPE III CLEAR	02/07/07		B. LLOYD	02/07/07
16oz (497.7cc ± 7.4cc OVERFLOW)				
FINISH: G.P.I. 48-400-405				
SCALE: FULL UNLESS NOTED				
DO NOT SCALE THIS DRAWING				
WHEATON®				
2001052114301198				
APPROX. WGT. 288.5 gr. 10.16oz				
AUTHORIZED C.J.M. 2/27/07				
REV	E.C.N. #	DATE		
B		02/07/07		
DWG. NO: C- W010514				