

-----NOTES-----

PRODUCT COMPONENTS

1. DIMENSIONS ARE IN INCHES AND [MILLIMETERS].
2. MILLIMETER DIMENSIONS ARE FOR REFERENCE ONLY.

ITEM	QTY	DESCRIPTION

TOLERANCE UNLESS NOTED: ±.010 [0.25]

6 PANELS
EQUALLY SPACED

∅1.616 [41.05]


LINER

1.051
[26.70]

CMA 38-430 FINISH
6 T.P.I. $\triangle B$

"H" .956-.938
[24.28-23.83]

THIRD ANGLE PROJECTION



∅1.638 [41.61]

"E" ∅1.404-1.383
[35.65-35.12]

"T" ∅1.491-1.477
[37.87-37.52]

THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION AND IS A COPYRIGHT OF DWK LIFE SCIENCES INC. IT MAY NOT BE DUPLICATED FOR ANY PURPOSE IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN CONSENT OF DWK LIFE SCIENCES INC. ALL RIGHTS RESERVED. ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF DWK LIFE SCIENCES SYSTEM ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

CAP MATERIAL: BLACK PHENOLIC
LINER MATERIAL: TAN PTFE FACED RUBBER $\triangle B C$

: \DWGS\240481

TITLE
38-430 CAP, PHEN BLK, PTFE/RUB LNR $\triangle B$

DATE DRAWN 03/31/09

PROJ#

DRAWN BY B. LLOYD

APPROX. WGT.

CAPACITY

SIZE

CHECKED BY TRJ 4/1/09

AUTHORIZED JB 4/1/09

FINISH CMA 38-430 $\triangle B$

DURAN WHEATON KIMBLE
1501 North 10th Street | Millville, NJ 08332-2038 | USA

-	CORPORATE NAME CHANGE	04/20/18 EL
C	2015090315124075	03/07/17 EL
B	2011101211223320	10/31/11 TP
A	2009033008031261	03/31/09 BL
REV	E.C.N. #	DATE

SCALE: FULL UNLESS NOTED. DO NOT SCALE THIS DRAWING.

DWG.NO:A- 240481