


:\00-089\W985872-REDWN-EXT
 THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION AND IS A COPYRIGHT OF WHEATON INDUSTRIES INC. IT MAY NOT BE DUPLICATED FOR ANY PURPOSE IN WHOLE OR IN PART WITHOUT THE EXPRESS WRITTEN CONSENT OF WHEATON INDUSTRIES INC. ALL RIGHTS RESERVED. ©

TOLERANCE: ±.005 [0.13] UNLESS NOTED


THIRD ANGLE PROJECTION




PRODUCT COMPONENTS

ITEM	QTY	DESCRIPTION

EXTERNAL SPECIFICATIONS



ACTUAL SIZE



SCALE: 2 X 1

UNCONTROLLED COPY

CRYULE MATERIAL: NATURAL POLYPROPYLENE
 CAP MATERIAL: NATURAL POLYPROPYLENE

STERILIZED

- 3. LOGO COLOR: MATTE WHITE.
- 2. APPROX WGT: CRYULE, 1.80 GRAMS; CAP, 0.51 GRAMS.
- 1. ALL DIMENSIONS ARE IN INCHES AND [MILLIMETERS].

TITLE 2ML EXT FS CRYOELITE™ NAT CAP PATCH STRL	
CAPACITY 3.2mL OVERFLOW	SIZE 2.0mL
FINISH 13mm SPECIAL	
SCALE: FULL UNLESS NOTED. DO NOT SCALE THIS DRAWING.	

DATE DRAWN 08/03/09	MAT'L
DRAWN BY B. LLOYD	APPROX. WGT.
CHECKED BY JB 8/12/09	AUTHORIZED 8/12/09
WHEATON Industries Inc. 1501 North 10th Street, Millville, NJ 08332-2038	

D	2012012615430034	04/10/12 BL
C	2009111310200232	01/07/10 BL
B	200909291442184	10/22/09 JB
A	20080409143405	08/03/09 BL
REV	E.C.N. #	DATE

DWG. NO: B- W985872-EXT