

PRODUCT SPECIFICATION

<i>Specification</i>	Product Code	W985917
	Description	CryoELITE® Cryogenic Vials
	Material	Polypropylene
	Colour	Red
	Autoclavable	No
	Pack quantity	500

<i>Dimension</i>	Height (mm)	50
	Diameter (mm)	12

DWK Life Sciences Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of DWK Life Sciences Limited by:

Signed: M.Ashley
Quality Technician
Date: 13 Aug 2020