

## PRODUCT SPECIFICATION

<i>Specification</i>	<b>Product Code</b>	W985921
	<b>Description</b>	CryoELITE® Cryogenic Vials
	<b>Material</b>	Polypropylene
	<b>Colour</b>	Blue
	<b>Autoclavable</b>	No
	<b>Pack quantity</b>	500

<i>Dimension</i>	<b>Height (mm)</b>	50
	<b>Diameter (mm)</b>	12

DWK Life Sciences Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of DWK Life Sciences Limited by:

Signed: M.Ashley  
Quality Technician  
Date: 13 Aug 2020