

Terumo Syringe

3-part Syringe (with or without Needle)



Product Solutions

Your syringes of choice

NEW



The Terumo logo consists of a red curved line above the word 'TERUMO' in a bold, green, sans-serif font.

Helping you
to help people




Offering
sustainable
solutions

Disposable syringes

for single use

The syringes are designed to be safe, convenient and represent the current state of the art in modern injection technology.

Carefully selected materials – free from PVC or natural rubber latex – together with a controlled production process result in products of excellent quality and performance.



Terumo Syringes and
Terumo Syringes with Needles

The range includes

- Concentric Luer Slip tip
- Eccentric Luer Slip tip
- Luer Lock tip



Smooth and perfectly controlled injection for reliable drug delivery

Made to support you in your daily practice

Every single detail is the result of Terumo's long time experience in this field which is appreciated by users throughout the world.

Performance enhancing features:

- siliconisation for optimal and smooth plunger movement
- firm and definite end-stop helps prevent the plunger from accidentally being pulled out
- finger grips support easy handling

Safety enhancing features:

- transparent barrel provides easy recognition of air bubbles
- clear printing allows precise dose management



ISO colour-coding of the needles:

- quick identification of the required size
- easy recognition of air bubbles or flashback through the translucent hub


Choose from various sizes optionally with needles

Terumo Syringes

Sterile 3-part syringes
in a choice of frequently
needed volumes

Terumo Syringes with needles

Sterile 3-part syringes
pre-connected with a sterile
hypodermic needle



Syringe in 1 ml volume

available with these needles:



25G Cannula (ø 0.5 mm)

Length 16 mm, Orange

26G Cannula (ø 0.45 mm)

Length 13 mm, Brown

Syringe in 2 ml volume

available with these needles:



21G Cannula (ø 0.8 mm)

Length 16/25/38 mm, Green

22G Cannula (ø 0.7mm)

Length 38 mm, Black

23G Cannula (ø 0.6 mm)


Length 25/32 mm, Blue

25G Cannula (ø 0.5 mm)

Length 16 mm, Orange

Syringe in 5 ml volume

available with these needles:



21G Cannula (ø 0.8 mm)

Length 38 mm, Green

22G Cannula (ø 0.7mm)


Length 38 mm, Black

23G Cannula (ø 0.6 mm)

Length 32 mm, Blue

Syringe in 10 ml volume

available with these needles:



20G Cannula (ø 0.9 mm)

Length 38 mm, Yellow

21G Cannula (ø 0.8 mm)

Length 38 mm, Green

22G Cannula (ø 0.7mm)

Length 38 mm, Black

For further ordering information, please refer to the Terumo Syringes /
Terumo Syringes with needle product information card.

Terumo Syringe

3-part Syringe (with or without Needle)



Dedicated to making things better

Terumo Syringes and Terumo Syringes with Needles will enhance Terumo's advanced Injection & Access portfolio

We understand our devices are an essential element for treating patients and saving lives – especially in crisis impacted scenarios. It is therefore our commitment to ensure stability of supply whilst bringing new products to the market.

With the new Terumo Syringes we provide solutions to medical institutions and caregivers.

Contributing to Society through Healthcare



Terumo Corporation
+81 3 6742 8500

Terumo Europe N.V.
+32 16 38 12 11

Your Terumo partner

www.terumo-europe.com

Please note that some products may not be available in all regions, check with your local Terumo contact.

© Registered Trademark
Published by Terumo Europe N.V.