


CATALOG NUMBER.....	SEE BELOW	
NOM. CAP. (uL).....	120	
MATERIAL.....	PS	
	INCH	MM
A (BASE LENGTH).....	5.03	127.8
B (BASE WIDTH).....	3.37	85.5
C (TOP LENGTH).....	4.87	123.7
D (TOP WIDTH).....	3.20	81.3
F (OVERALL HEIGHT).....	0.57	14.4
G DIA.....	0.15	3.7
H (A-1 LOCATION) (X).....	0.48	12.1
J (A-1 LOCATION) (Y).....	0.35	9.0
K (WELL-TO-WELL SPACING).....	.18	4.5
L (WELL DEPTH).....	0.46	11.7
M (BOTTOM OF FILM).....	.07	1.7
*P (LID LENGTH).....	5.02	127.4
*Q (LID WIDTH).....	3.36	85.3
*R (LID HEIGHT).....	.36	9.1
S (STACKED HEIGHT - NO LID).....	1.09	27.7
*T (STACKED HEIGHT W/LIDS).....	1.33	33.8
*U (ASSEMBLY HEIGHT W/LID).....	.665	16.9

CATALOG NO.	SURFACE	COLOR	STERILE	WITH LID*	UNITS/PK/CASE
142761	CELL CULTURE	BLACK	YES	YES	10/30
142762	CELL CULTURE	WHITE	YES	YES	10/30
164730	CELL CULTURE	BLACK	YES	NO	10/30
152041	COLLAGEN I	BLACK	NO	YES	5/20
152029	POLY-D-LYSINE	BLACK	NO	YES	5/20
242763	NON-TREATED	WHITE	NO	NO	10/30
242764	NON-TREATED	BLACK	NO	NO	10/30

* DIM.S APPLY ONLY WHERE "WITH LID" = YES, USING NUNC LID #8337.

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher
SCIENTIFIC**

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

DATE	JDD	APP'D	PVC	PSS 2427XX	REV.	0207
CHK'D	PVC	2/16/07				

DESCRIPTION: 384 WELL, PS, OPTICAL BOTTOM PLATE