


CATALOG NUMBER	SEE BELOW	
WORKING RANGE (µL)	20 – 450µL	
MATERIAL	PP	
	MM	INCH
A. BASE LENGTH	127,8	5,03
B. BASE WIDTH	85,7	3,37
C. TOP LENGTH	123	4,84
D. TOP WIDTH	80,9	3,1
E. FLANGE LONG SIDE	7,5	0,30
F. OVERALL HEIGHT	14,5	0,58
G1, DIA TOP	Ø8,5	Ø0,32
G2, DIA BOTTOM	Ø7,7	Ø0,30
H. A - I LOCATION (X)	14,3	0,56
J. A - I LOCATION (Y)	11,3	0,45
K. WELL TO WELL SPACING	9	0,35
L. WELL DEPTH	11,2	0,44
M. BOTTOM OF WELL DISTANCE	3,1	0,12
N. FLANGE SHORT SIDE	2,5	0,10

CATALOG NO.	SURFACE	COLOR	STERILE	WITH LID*	UNITS/ PKG/CASE
267245	NON-TREATED	NATURAL	-	-	10/120
267334	NON-TREATED	NATURAL	-	-	10/120
267342	NON-TREATED	BLACK	-	-	10/120
267350	NON-TREATED	WHITE	-	-	10/120
267369	NON-TREATED	RED	-	-	10/120
267385	NON-TREATED	BLUE	-	-	10/120
267407	NON-TREATED	YELLOW	-	-	10/120

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
  - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

## PLATE



**ThermoFisher**  
SCIENTIFIC

Roskilde Site • Kamstrupvej 90 • P.O.Box 280  
DK-4000 Roskilde • Denmark  
Tel +45 4631 2000 • Fax +45 4631 2042  
E-mail: info.nunc@thermofisher.com

[www.nuncbrand.com](http://www.nuncbrand.com)

SUBJECT: U96 MICROWELL PP.		DRAWN BY/DATE:	CHANGE ORDER NO.:	VALID FROM:	APPROVED BY/DATE:	FORMAT:
ITEM NUMBER:	MATERIAL:	DRAWING NO.:	VERSION:	PAGE TYPE:	SCALE:	METRIC:

CONFIDENTIAL  
THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESS-PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC