



Roskilde Site • Kamstrupvej 90 • P.O.Box 280
 DK-4000 Roskilde • Denmark
 Tel +45 4631 2000 • Fax +45 4631 2042
 E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

SUBJECT: Stopper for CryoTube. Ersattes af 2308-C.idw	DRAWN BY/DATE: AVG/2007.07.14	CHANGE ORDER NO.: CO-3989	VALID FROM: 2011.07.07	APPROVED BY/DATE: ROSPSI 2011.07.04	FORMAT: A4	
	ITEM NUMBER: -	MATERIAL: PP	DRAWING NO.: 2308	VERSION: 12	PAGE TYPE: -C	SCALE: 5:1