


- Useful in all areas of cell culture including scale-up and cloning
- Raised well rims to lower risk of cross contamination
- Excellent optical quality
- Nunclon<sup>™</sup>Δ certified
- Multidish 6 is available with Poly-D-Lysine and Collagen Coated surface. Produced in clean environment

**Multidishes**

Polystyrene. With lid




Cat. No.	176740*	176597	140675**	140685**	152035	152034	176600	150628
Number of wells	4	4	6	6	6	6	8	12
Surface	Nunclon Δ	Nunclon Δ	Nunclon Δ	Nunclon Δ	Poly-D-Lysine	Collagen 1, rat tail	Nunclon Δ	Nunclon Δ
Bottom shape	Flat	Flat	Flat	Flat	Flat	Flat	Flat	Flat
Sterile	+	+	+	+	-	-	+	+
Suggested working vol., ml/well	1	5	3	3	3	3	3	2
Culture area, cm <sup>2</sup> /well	1.9	16.0	9.6	9.6	9.6	9.6	8.6	3.5
External dimensions, mm	66 x 66	132 x 32	128 x 86	128 x 86	128 x 86	128 x 86	132 x 82	128 x 86
Units per pack/case	4/120	10/100	1/75	5/85	5/20	5/20	10/100	1/75

As from 2003 \*No. 176740 has a New Design. \*\*New MD 6: No. 140675 and 140685 replace No. 152795 and 150229.

**Multidishes Nunclon<sup>™</sup>Δ**

Polystyrene. With lid. Sterile



Cat. No.	142475**	142485**	150687	152640
Number of wells	24	24	48	48
Bottom shape	Flat	Flat	Flat	Flat
Suggested working vol., ml/well	1.0	1.0	0.5	0.5
Culture area, cm <sup>2</sup> /well	1.9	1.9	1.1	1.1
External dimensions, mm	128 x 86	128 x 86	128 x 86	128 x 86
Units per pack/case	1/75	5/85	1/75	5/85

\*\*New MD 24: No. 142475 and 142485 replace No. 143982 and 146485

**LITERATURE:**

Nunc Bulletin No: 2, 3

Page 293

**ACCESSORIES:**

Cell Scraper

Page 15

Serological Pipettes

Page 15

Nunc CC Inserts

Page 24