



CATALOG NUMBER.....	376813	
NOM. CAP. (mL).....	200	
CLOSURE SIZE (mm).....	28	
THREAD (T.P.I.).....	6	
TUBE.....	PP	
CLOSURE.....	HOPE	
A (NECK I.O.).....	INCH	MM
B (HT. W/ CLOSURE).....	1.02	25.9
C (HT. W/O CLOSURE).....	5.40	137.2
D (O.I.A.).....	5.38	136.7
(W.T.) WALL THICKNESS:	2.36	59.9
	.04	1.0

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Ponderosa Creek Drive
Rochester, NY 14625-2385, USA

NAME	MC ²	MP'D	PVC	REV.	0102
FOR	1/15/02	1/15/02			
DESCRIPTION:	PSS 376813				

200mL CONICAL TUBE, PP

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.