



Gravering af tal fra 1 og fortløbende, 0.03mm gravureddybde, talhøjde 1.2mm, bemærk en bundbrik skal ikke graves.

Tolerances unless otherwise specified:
 Without decimal: ± 0.5 mm
 1 place decimal: ± 0.1 mm
 2 place decimal: ± 0.05 mm

SUBJECT:

Prøverør $\varnothing 13 \times 100$

SCALE:

2:1(1:1)

COMPLETE BY/DATE:

Tbs/15.09.99

VERSION:

2.0

APPROVED BY/DATE:

T / 15.09.99

MATERIAL:

PS

VALID FROM:

15.09.99

DRAWING:

1408

PAGE

1

OF

1



Nunc A/S
 Kamstrupvej 90
 P.O. Box 280
 Dk-4000 Roskilde
 Denmark

Tel 45 46 35 90 65
 Fax 45 46 32 06 99
 e-mail: infacety@nunc.dk

A Subsidiary of
 Sybron
 International

Nunc Nunc
 International