




CATALOG NO.....	178599	
	INCH	MM
A (LENGTH).....	2.98	75.7
B (WIDTH).....	1.02	25.9
C (O/A HEIGHT).....	.51	12.9
D.....	.22	5.6
E.....	.23	5.8
F.....	.356	9.05
G.....	.356	9.05

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Paragona Creek Drive
Rochester, NY 14625-2383, USA

DATE: JD P.V.C.
8/25/08 8/25/08
REV: 0808

PSS - 178599

16-WELL CHAMBER SLIDE