


2:1



Tolerances unless otherwise specified:
 Without decimal: ± 0.5 mm
 1 place decimal: ± 0.1 mm
 2 place decimal: ± 0.05 mm

Undefined radius: R 0.20 mm

ThermoFisher
 S C I E N T I F I C

SUBJECT:
 Connector for Cell factory
 Cat. No. 73330

Change Order no.:	00846	Changed by/date:	TBS 2004.09.09	Approved by/date:	PS 2004.09.09	Valid from:	2004.09.09
Change Order no.:	00087	Changed by/date:	TBS/R 2004.07.26	Approved by/date:	PS 2004.07.26	Valid from:	2004.07.26

www.nuncbrand.com
 Roskilde Site - Kamstrupvej 90 • P.O.Box 280
 DK-4000 Roskilde • Denmark
 Tel +45 4631 2000 • Fax +45 4631 2042
 E-mail: info.nunc@thermofisher.com

DRAWN BY/DATE:	JR/MS/98.11.24	CHANGE ORDER NO.:	CO-1181	VALID FROM:	2009.06.17	APPROVED BY/DATE:	JR2009.06.17	FORMAT:	A3
ITEM NUMBER:	73330	MATERIAL:	White Teflon	DRAWING NO.:	3033	VERSION:	10	PAGE TYPE:	-M
							SCALE:	5:1 (2:1)	METRIC