



A



A

A-A



2:1

Tolerances unless otherwise specified:
 Without decimal: ± 0.5 mm
 1 place decimal: ± 0.1 mm
 2 place decimal: ± 0.05 mm

Preliminary drawing

SUBJECT:

Module F16 mounted in frame.

Resille Site - Kamintrasse 90 - P.O.Box 280
 DK-4000 Roskilde - Denmark
 Tel: +45 4631 2000 - Fax: +45 4631 2042
 E-mail: info.nunc@thermofisher.com

ThermoFisher
 SCIENTIFIC

www.nuncbrand.com

SCALE:	1:1	COMPLETED BY/DATE:	PS/95.10.01
VERSION:	1.0	APPROVED BY/DATE:	PS/95.10.01
MATERIAL:	PS	VALID FROM:	95.10.01
DRAWING NO.:	2132	PAGE	OF
			1 1