



Section CC

All dimensions are strictly related to the combination of strips mounted in the frame.

Tolerances unless otherwise specified:
 Without decimal: ± 0.5 mm
 1 place decimal: ± 0.1 mm
 2 place decimal: ± 0.05 mm

ThermoFisher
 S C I E N T I F I C

Roskilde Site • Kamstrupvej 90 • P.O.Box 280
 DK-4000 Roskilde • Denmark
 Tel +45 4631 2000 • Fax +45 4631 2042
 E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

SUBJECT: LockWell C8 Module. Assembly drawing.		DRAWN BY/DATE: BGL/25.08.94		CHANGE ORDER NO.:	VALID FROM:	APPROVED BY/DATE:	FORMAT:
		ITEM NUMBER:	MATERIAL:	CO-2089	2007.12.07	JR 2007.12.04	A3
			PS/ABS	DRAWING NO.:	VERSION:	PAGE TYPE:	METRIC
				2459	11	-C	Scale: 1:1, 2:1