



ThermoFisher
SCIENTIFIC

Roskilde Site - Kamstrupvej 90 - P.O.Box 280
DK-4000 Roskilde - Denmark
Tel +45 4631 2000 • Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

Ⓐ & Ⓑ to be measured 4.2 mm from the edge. At the chamfered corner 10.5 mm from the edge.

Tolerances unless otherwise specified:
Without decimal: ±0.5 mm
1 place decimal: ±0.1 mm
2 place decimal: ±0.05 mm

SUBJECT: Frame for LockWell		ITEM NUMBER: MS/JR/13.06.02	MATERIAL: ABS	CHANGE ORDER NO.:	VERSION: 10	APPROVED BY/DATE: JR 2007.11.12	FORMAT: A3
		VALID FROM: 2007.11.19	DRAWING NO.:	SCALE: 1:1	PAGE TYPE: -C	METRIC ASSOCIATION	
			2370				