



1:1 (Lid + plate)

Tolerances unless otherwise specified:
 Without decimal: ±0.5 mm
 1 place decimal: ±0.1 mm
 2 place decimal: ±0.05 mm
 Undefined radius: R 0.20 mm

ThermoFisher
 S C I E N T I F I C

Nunc AS - Karstrupvej 90 - P.O. Box 280
 DK-4000 Roskilde - Denmark
 Tel +45 4631 2000 - Fax +45 4631 2042
 E-mail: info.nunc@thermofisher.com

www.thermoscientific.com

SUBJECT: Lid for MiniTray plate
 (Lid + MiniTray Plate)

DRAWING/ORDER NO.:	CHANGE ORDER NO.:	VALID FROM:	APPROVED BY/DATE:	FORMAT:
98E2012-03-15	00-4446	2012.03.20	ROSKGA 2012.03.15	A3
ITEM NUMBER:	MATERIAL:	DRAWING NO.:	PAGE TYPE:	METRIC
PS		2145	-C	SCALE: 2x(1:5)