


CATALOG NUMBER	SEE BELOW	
WORKING RANGE (µL)	50 – 250 µL	
MATERIAL	PS	
	MM	INCH
A, BASE LENGTH	127,7	5
B, BASE WIDTH	85,6	3,38
C, TOP LENGTH	123,3	4,84
D, TOP WIDTH	81,2	3,20
E, FLANGE LONG SIDE	2,7	0,11
F, OVERALL HEIGHT	14,6	0,57
G1, DIA. TOP	Ø7,05	Ø0,28
G2, DIA. BOTTOM	Ø6,55	Ø0,26
H, A - I LOCATION (X)	14,4	0,55
J, A - I LOCATION (Y)	14,4	0,45
K, WELL TO WELL SPACING	9	0,35
L, WELL DEPTH	11,2	0,44
M, BOTTOM OF WELL DISTANCE	3,6	0,14
N, FLANGE SHORT SIDE	2,2	0,09
*P, LID LENGTH	127,4	5,02
*Q, LID WIDTH	85,4	3,36
*R, LID HEIGHT	9	0,35
S, STACKED HEIGHT – NO LID	27	1,06
*T, STACKED HEIGHT W/LIDS	32,3	1,27
*U, ASSEMBLY HEIGHT W/LID	15,7	0,62

CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	UNITS/PK/CASE	CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	UNITS/PK/CASE
136101	CELL CULTURED	WHITE	+	+	1/50	237108	NON-TREATED	BLACK	-	-	10/180
136102	CELL CULTURED	WHITE	+	+	10/160	436110	MAXISORP	WHITE	-	-	10/80
137101	CELL CULTURED	BLACK	+	+	1/50	437111	MAXISORP	BLACK	-	-	10/80
137103	CELL CULTURED	BLACK	+	+	10/160	436111	POLYSORP	WHITE	-	-	10/80
236105	NON-TREATED	WHITE	+	+	1/50	437112	POLYSORP	BLACK	-	-	10/80
236107	NON-TREATED	WHITE	+	+	10/160						
236108	NON-TREATED	WHITE	-	-	10/180						
237105	NON-TREATED	BLACK	+	+	1/50						
237107	NON-TREATED	BLACK	+	+	10/160						

\* DIM.S APPLY ONLY WHERE "WITH LID" = +. USING NUNC LID 2205

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
  - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



<b>ThermoFisher</b> <b>SCIENTIFIC</b>		Roskilde Site • Kamstrupvej 90 • P.O.Box 280 DK-4000 Roskilde • Denmark Tel +45 4631 2000 • Fax +45 4631 2042 E-mail: Info.nunc@thermofisher.com		<a href="http://www.nuncbrand.com">www.nuncbrand.com</a>	
SUBJECT:	DRAWN BY/DATE:	CHANGE ORDER NO.:	VALID FROM:	APPROVED BY/DATE:	FORMAT:
F96 MICROWELL PLATE, PS					
ITEM NUMBER:	MATERIAL:	DRAWING NO.:	VERSION:	PAGE TYPE:	SCALE:
					METRIC

CONFIDENTIAL  
 THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESS-PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC