


CATALOG NUMBER	SEE BELOW	
WORKING RANGE ( $\mu$ L)	50 – 250 $\mu$ L	
MATERIAL	PS	
	MM	INCH
A. BASE LENGTH	127,7	5,03
B. BASE WIDTH	85,7	3,37
C. TOP LENGTH	123,5	4,86
D. TOP WIDTH	81,3	3,20
E. FLANGE LONG SIDE	2,7	0,11
F. OVERALL HEIGHT	14,5	0,57
G1. DIA TOP	$\varnothing$ 6,6	$\varnothing$ 0,26
G2. DIA BOTTOM	$\varnothing$ 6,1	$\varnothing$ 0,24
H. A - I LOCATION (X)	14,3	0,56
J. A - I LOCATION (Y)	11,4	0,45
K. WELL TO WELL SPACING	9	0,35
L. WELL DEPTH	11	0,43
M. BOTTOM OF WELL DISTANCE	3,5	0,14
N. FLANGE SHORT SIDE	2,7	0,11

CATALOG NO	SURFACE	COLOR	STERILE	WITH LID	UNITS/PKCASE
437796	MAXISORP	WHITE	-	-	5/60
437842	POLYSORP	WHITE	-	-	5/60

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
  - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



**Thermo Fisher**  
SCIENTIFIC

Roskilde Site • Kamstrupvej 90 • P.O. Box 280  
DK-4000 Roskilde • Denmark  
Tel +45 4631 2000 • Fax +45 4631 2042  
E-mail: info.nunc@thermofisher.com

[www.nuncbrand.com](http://www.nuncbrand.com)

SUBJECT:	DRAWN BY/DATE:	CHANGE ORDER NO.:	VALID FROM:	APPROVED BY/DATE:	FORMA:
C96 MICROWELL PLATE, WHITE, PS					
ITEM NUMBER:	MATERIAL:	DRAWING NO.:	VERSION:	PAGE TYPE:	SCALE:
					METRIC

CONFIDENTIAL  
THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC