



© 201




A

10:1



A-A



Tolerances unless otherwise specified:  
Without decimal:  $\pm 0.5$  mm  
1 place decimal:  $\pm 0.1$  mm  
2 place decimal:  $\pm 0.05$  mm  
Undefined radius: R 0.20 mm

**ThermoFisher**  
SCIENTIFIC

Roskilde Site - Kamstrupvej 90 - P.O. Box 280  
DK-4000 Roskilde - Denmark  
Tel +45 4631 2000 - Fax +45 4631 2042  
E-mail: info.nunc@thermofisher.com  
[www.nuncbrand.com](http://www.nuncbrand.com)

SUBJECT: Loop for 10 micro liter		DRAWN BY/DATE: TBS/24.02.00	CHANGE ORDER NO.: CO-1181	VALID FROM: 2008.05.08	APPROVED BY/DATE: HES 2008.05.05	FORMAT: A3
		ITEM NUMBER:	DRAWING NO.: 2158	VERSION: 10	PAGE TYPE: -C	SCALE: 2:1(10:1)
		MATERIAL: PS	METRIC UNITS			