


| CATALOG NUMBER | SEE BELOW | |
|----------------------------|------------|-------|
| WORKING RANGE (µL) | 50 - 250µL | |
| MATERIAL | PS | |
| | MM | INCH |
| A. BASE LENGTH | 127,8 | 5,03 |
| B. BASE WIDTH | 85,5 | 3,37 |
| C. TOP LENGTH | 123,3 | 4,84 |
| D. TOP WIDTH | 81,2 | 3,20 |
| E. FLANGE LONG SIDE | 7,5 | 0,3 |
| F. OVERALL HEIGHT | 14,4 | 0,57 |
| G1. DIA TOP | Ø7 | Ø0,28 |
| G2. DIA BOTTOM | Ø6,2 | Ø0,24 |
| H. A - I LOCATION (X) | 14,3 | 0,56 |
| J. A - I LOCATION (Y) | 11,3 | 0,45 |
| K. WELL TO WELL SPACING | 9 | 0,35 |
| L. WELL DEPTH | 11,4 | 0,45 |
| M. BOTTOM OF WELL DISTANCE | 3,2 | 0,13 |
| N. FLANGE SHORT SIDE | 2,4 | 0,1 |
| *P. LID LENGTH | 127,4 | 5,02 |
| *Q. LID WIDTH | 85,4 | 3,36 |
| *R. LID HEIGHT | 9 | 3,36 |
| S. STACKED HEIGHT - NO LID | 26,9 | 1,06 |
| *T. STACKED HEIGHT W/LIDS | 31,2 | 1,23 |
| *U. ASSEMBLY HEIGHT W/LID | 15,6 | 0,61 |

| CATALOG NO | SURFACE | COLOR | STERILE | WITH LID* | BAR CODE | UNITS/ PK/CASE | CATALOG NO | SURFACE | COLOR | STERILE | WITH LID* | BAR CODE | UNITS/ PK/CASE |
|------------|---------------|-------|---------|-----------|----------|----------------|------------|-------------|-------|---------|-----------|----------|----------------|
| 156545 | CELL CULTURE | CLEAR | + | - | - | 10/180 | 260844 | NON-TREATED | CLEAR | - | - | - | 10/180 |
| 161093 | CELL CULTURE | CLEAR | + | - | - | 10/180 | 260860 | NON-TREATED | CLEAR | + | + | - | 10/160 |
| 167008 | CELL CULTURE | CLEAR | + | - | - | 1/50 | 260887 | NON-TREATED | CLEAR | + | + | + | 10/160 |
| 168025 | CELL CULTURE | CLEAR | + | - | - | 1/50 | 260895 | NON-TREATED | CLEAR | - | + | - | 10/160 |
| 152038 | COLLAGEN I | CLEAR | - | - | - | 5/20 | 456529 | POLYSORP | CLEAR | - | - | - | 10/180 |
| 152039 | POLY-D-lysine | CLEAR | - | - | - | 5/20 | 456537 | MAXISORP | CLEAR | - | - | - | 10/180 |
| 178747 | CC™ | CLEAR | + | - | - | 5/50 | 460984 | MAXISORP | CLEAR | - | - | + | 10/180 |
| 258510 | NON-TREATED | CLEAR | + | + | - | 10/180 | 467320 | MEDISORP | CLEAR | - | - | - | 5/60 |
| 260836 | NON-TREATED | CLEAR | - | - | - | 10/180 | 467340 | MULTISORP | CLEAR | - | - | - | 5/60 |

* DIM.S APPLY ONLY WHERE "WITH LID" = +. USING NUNC LID 2205

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
 - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



Roskilde Site • Kamstrupvej 90 • P.O.Box 280
DK-4000 Roskilde • Denmark
Tel +45 4631 2000 • Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

| | | | | | | |
|---|--|----------------|-------------------|--------------|-------------------|-------------------------------|
| SUBJECT: | | DRAWN BY/DATE: | CHANGE ORDER NO.: | VALID FROM: | APPROVED BY/DATE: | FORM |
| F96 MICROWELL PLATE, HIGH FLANGE, PS | | ITEM NUMBER: | MATERIAL: | DRAWING NO.: | VERSION: | PAGE TYPE: SCALE: METRIC/INCH |

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC