


CATALOG NUMBER	SEE BELOW	
WORKING RANGE (µL)	50 - 250µL	
MATERIAL	PS	
	MM	INCH
A. BASE LENGTH	127,8	5,03
B. BASE WIDTH	85,5	3,37
C. TOP LENGTH	123,3	4,84
D. TOP WIDTH	81,2	3,20
E. FLANGE LONG SIDE	7,5	0,3
F. OVERALL HEIGHT	14,4	0,57
G1. DIA TOP	Ø7	Ø0,28
G2. DIA BOTTOM	Ø6,2	Ø0,24
H. A - I LOCATION (X)	14,3	0,56
J. A - I LOCATION (Y)	11,3	0,45
K. WELL TO WELL SPACING	9	0,35
L. WELL DEPTH	11,4	0,45
M. BOTTOM OF WELL DISTANCE	3,2	0,13
N. FLANGE SHORT SIDE	2,4	0,1
*P. LID LENGTH	127,4	5,02
*Q. LID WIDTH	85,4	3,36
*R. LID HEIGHT	9	3,36
S. STACKED HEIGHT - NO LID	26,9	1,06
*T. STACKED HEIGHT W/LIDS	31,2	1,23
*U. ASSEMBLY HEIGHT W/LID	15,6	0,61

CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	BAR CODE	UNITS/ PK/CASE	CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	BAR CODE	UNITS/ PK/CASE
156545	CELL CULTURE	CLEAR	+	-	-	10/180	260844	NON-TREATED	CLEAR	-	-	-	10/180
161093	CELL CULTURE	CLEAR	+	-	-	10/180	260860	NON-TREATED	CLEAR	+	+	-	10/160
167008	CELL CULTURE	CLEAR	+	-	-	1/50	260887	NON-TREATED	CLEAR	+	+	+	10/160
168025	CELL CULTURE	CLEAR	+	-	-	1/50	260895	NON-TREATED	CLEAR	-	+	-	10/160
152038	COLLAGEN I	CLEAR	-	-	-	5/20	456529	POLYSORP	CLEAR	-	-	-	10/180
152039	POLY-D-lysine	CLEAR	-	-	-	5/20	456537	MAXISORP	CLEAR	-	-	-	10/180
178747	CC™	CLEAR	+	-	-	5/50	460984	MAXISORP	CLEAR	-	-	+	10/180
258510	NON-TREATED	CLEAR	+	+	-	10/180	467320	MEDISORP	CLEAR	-	-	-	5/60
260836	NON-TREATED	CLEAR	-	-	-	10/180	467340	MULTISORP	CLEAR	-	-	-	5/60

\* DIM.S APPLY ONLY WHERE "WITH LID" = +. USING NUNC LID 2205

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
  - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



Roskilde Site • Kamstrupvej 90 • P.O.Box 280  
DK-4000 Roskilde • Denmark  
Tel +45 4631 2000 • Fax +45 4631 2042  
E-mail: info.nunc@thermofisher.com

[www.nuncbrand.com](http://www.nuncbrand.com)

SUBJECT:		DRAWN BY/DATE:	CHANGE ORDER NO.:	VALID FROM:	APPROVED BY/DATE:	FORM
F96 MICROWELL PLATE, HIGH FLANGE, PS		ITEM NUMBER:	MATERIAL:	DRAWING NO.:	VERSION:	PAGE TYPE:
		SCALE:	METRIC	UNITS	SCALE	SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC