


CATALOG NUMBER.....	SEE BELOW	
NOM. CAP. (uL).....	120	
MATERIAL.....	PS	
	INCH	MM
A (BASE LENGTH).....	5.03	127.8
B (BASE WIDTH).....	3.37	85.5
C (TOP LENGTH).....	4.87	123.7
D (TOP WIDTH).....	3.20	81.3
E (FLANGE LONG SIDE).....	0.25	6.4
F (OVERALL HEIGHT).....	0.57	14.4
G DIA.....	0.15	3.7
H (A-1 LOCATION) (X).....	0.48	12.1
J (A-1 LOCATION) (Y).....	0.35	9.0
K (WELL-TO-WELL SPACING).....	0.18	4.5
L (WELL DEPTH).....	0.46	11.7
M (BOTTOM OF WELL DISTANCE).....	.07	1.7
N (FLANGE SHORT SIDE).....	.10	2.5
*P (LID LENGTH).....	5.02	127.4
*Q (LID WIDTH).....	3.36	85.3
*R (LID HEIGHT).....	.36	9.1
S (STACKED HEIGHT - NO LID).....	1.09	27.7
*T (STACKED HEIGHT W/LIDS).....	1.33	33.8
*U (ASSEMBLY HEIGHT W/LID).....	.665	16.9

CATALOG NO.	SURFACE	COLOR	STERILE	WITH LID*	UNITS/PK/CASE	CATALOG NO.	SURFACE	COLOR	STERILE	WITH LID*	UNITS/PK/CASE
164564	CELL CULTURE	BLACK	YES	YES	10/30	242765	NON-TREATED	CLEAR	YES	YES	1/30
164610	CELL CULTURE	WHITE	YES	YES	10/30	262260	NON-TREATED	BLACK	NO	NO	25/100
164688	CELL CULTURE	CLEAR	YES	YES	10/30	262360	NON-TREATED	WHITE	NO	NO	25/100
164571	CELL CULTURE	BLACK	YES	NO	25/100	262160	NON-TREATED	CLEAR	NO	NO	25/100
165195	CELL CULTURE	WHITE	YES	NO	25/100	265202	NON-TREATED	CLEAR	YES	YES	25/100
164555	CELL CULTURE	CLEAR	YES	NO	25/100	265203	NON-TREATED	CLEAR	YES	NO	25/100
264556	NON-TREATED	BLACK	NO	NO	10/30	460435	POLYSDRPP®	BLACK	NO	NO	10/30
264572	NON-TREATED	WHITE	NO	NO	10/30	460518	MAXISDRP™	BLACK	NO	NO	10/30
265196	NON-TREATED	CLEAR	NO	NO	10/30	460372	MAXISDRP™	WHITE	NO	NO	10/30
242757	NON-TREATED	CLEAR	YES	YES	10/30	464718	MAXISDRP™	CLEAR	NO	NO	10/30

\* DIM.S APPLY ONLY WHERE "WITH LID" = YES, USING NUNC LID #8337.

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

DATE: J00	APP'D: GRP	PSS 1646XX	REV: 0706
CHK'D: AEG	7/24/06		
DESCRIPTION:			

384 WELL, PS, FLAT BOTTOM PLATE