





| | | |
|---------------------------------|-----------|--------|
| CATALOG NUMBER..... | SEE BELOW | |
| NOM. CAP. (uL)..... | 400 | |
| MATERIAL..... | PS | |
| | INCH | MM |
| A (BASE LENGTH)..... | 5.030 | 127.76 |
| B (BASE WIDTH)..... | 3.365 | 85.47 |
| C (TOP LENGTH)..... | 4.87 | 123.3 |
| D (TOP WIDTH)..... | 3.20 | 81.8 |
| E (FLANGE HEIGHT)..... | 0.25 | 6.2 |
| F (OVERALL HEIGHT)..... | 0.59 | 14.9 |
| G1 DIA. (BOTTOM OF WELL)..... | .249 | 6.32 |
| G2 DIA. (TOP OF WELL)..... | .254 | 6.45 |
| H (A-1 LOCATION) (X)..... | 0.57 | 14.4 |
| J (A-1 LOCATION) (Y)..... | 0.44 | 11.2 |
| K (WELL-TO-WELL SPACING)..... | 0.35 | 9.0 |
| L (WELL DEPTH)..... | 0.48 | 12.2 |
| M (BOTTOM OF GLASS)..... | .09 | 2.3 |
| *N (LID LENGTH)..... | 5.02 | 127.4 |
| *P (LID WIDTH)..... | 3.36 | 85.3 |
| *Q (LID HEIGHT)..... | .36 | 9.0 |
| *R (STACKED HEIGHT W/LIDS)..... | 1.31 | 33.3 |
| *S (ASSEMBLY HEIGHT W/LID)..... | .65 | 16.6 |

| CATALOG NUMBER | SURFACE | COLOR | STERILE | WITH LID* | UNITS PER PACK/CASE |
|----------------|--------------|-------|---------|-----------|---------------------|
| 164588 | CELL CULTURE | BLACK | YES | YES | 6/30 |
| 164590 | CELL CULTURE | WHITE | YES | YES | 6/30 |
| 160376 | CC* | BLACK | YES | YES | 6/30 |
| 160377 | CC* | WHITE | YES | YES | 6/30 |
| 265300 | NON-TREATED | BLACK | NO | NO | 5/30 |

*DIM'S APPLY ONLY WHERE "WITH LID" = YES. USING NUNC LID #9326B

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panoroma Creek Drive
Rochester, NY 14625-2385, USA

| | | | | | |
|------|-----|-------|--------|------------|------|
| DATE | JDD | APP'D | PVC | | REV. |
| DATE | KDW | DATE | 6/8/06 | PSS 1645XX | 0506 |

DESCRIPTION: 96 WELL, PS, GLASS BOTTOM PLATE