


|                                      |            |       |
|--------------------------------------|------------|-------|
| CATALOG NUMBER.....                  | 153732     |       |
| NOM. CAP (mL).....                   | 260        |       |
| CAP TYPE.....                        | VENT/CLOSE |       |
| CULTURE AREA (cm <sup>2</sup> )..... | 80         |       |
| MATERIAL                             | PS         |       |
| FLASK                                | BLUE HDPE  |       |
| CAP                                  |            |       |
|                                      |            |       |
| A (NECK ID.).....                    | INCH       | MM    |
| B (LENGTH W/O CAP).....              | 0.77       | 19.7  |
| C (LENGTH W/ CAP).....               | 5.62       | 142.7 |
| D (BODY WIDTH).....                  | 5.72       | 145.3 |
| E (BODY DEPTH).....                  | 3.08       | 78.2  |
|                                      | 1.44       | 36.6  |



SPECS ON BACK

NOT TO SCALE

**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

|      |     |      |         |     |      |
|------|-----|------|---------|-----|------|
| FASE | MC2 | DATE | PVC     | REV | 0505 |
| CRD  | RDR |      | 6/29/05 |     |      |

PSS 153732

80cm<sup>2</sup> NUNC LON FLASK, STERILE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.