


CATALOG NUMBER.....:		156502	
NOM. CAP. (mL).....:		800	
NOM. TYPE.....:		VENT/CLOSE	
CULTURE AREA (cm ²).....:		175	
MATERIAL	FLASK.....:	PS	
	CAP.....:	BLUE HDPE	
		INCH	MM
A (NECK I.D.).....:		0.99	25.1
B (LENGTH W/O CAP).....:		7.64	194.1
C (LENGTH W/ CAP).....:		7.59	192.8
D (BODY WIDTH).....:		4.75	120.7
E (BODY DEPTH).....:		2.01	51.1



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC ²	APP'D	PVC	PSS 156502	REV.	0505
CK'D	RDR	6/29/05				
DESCRIPTION: 175cm ² NUNCLON FLASK, STERILE						