



ThermoFisher S C I E N T I F I C MD24 2nd Generation		ThermoFisher Scientific Esbjergvej 60 • 8500 St. Thors DK-8500 Roskilde, Denmark Tel. +45 4671 2000 • Fax. +45 4671 7042 Email: thermo@thermo.com	www.thermo.com www.nuncbrand.com
MODEL NUMBER MD24-020-1	PART NUMBER 2097424-1	PRODUCT CODE AN059103413	UNIT OF MEASURE 10 ea
MATERIAL PS	FINISH 10 ea	WEIGHT 3235	DIMENSIONS L x W x H 111 x 111 x 32

Remarks:
 Ⓣ Measured in the range 0 to 12,70 mm. from theoretical corner.
 Ⓣ Overall tolerance: ± 0.05 mm.
 Ⓣ Measured from the center of the well. Tolerance includes warp and twist of well and part.