


Text font: Arial
Text size: 3 mm
0,2 positive


Text font: Arial
Text size: 2 mm
0,2 positive

NUNC logo to be polished
Text font: Arial
Text size: 3 mm
0,05 positive

A (4 : 1)



Cavity no.
Text font: Arial
Text size: 1mm
0,05 positive

 = Surfacefrosting, charmilles 25

Undefined radius: R 0.2 mm
Tolerances unless otherwise specified:
Without decimal: ±0.5 mm
1 place decimal: ±0.1 mm
2 place decimal: ±0.05 mm

Thermo Fisher
S C I E N T I F I C

Nunc A/S · Kamstrupvej 90 · P.O.Box 280
DK-4000 Roskilde · Denmark
Tel +45 4631 2000 · Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.thermoscientific.com

SUBJECT: MD24 2nd generation		DRAWN BY/DATE: AVG/12.03.07	CHANGE ORDER NO.: CO-2859	VALID FROM: 2011-06-16	APPROVED BY/DATE: ROSPSI 2011-05-30	FORMAT: A3
ITEM NUMBER:	MATERIAL: PS	DRAWING NO.: 3225	VERSION: 10-a	PAGE TYPE: -C	SCALE: 1:1	