



Tolerances unless otherwise specified:

Without decimal: $\pm 0.5\text{mm}$

1 place decimal: $\pm 0.1\text{mm}$

2 place decimal: $\pm 0.05\text{mm}$

Change Order no.:	00641	Changed by/date:	TBS 25.03.04	Approved by/date:	PS /25.03.04	Valid from:	25.03.04
-------------------	-------	------------------	-----------------	-------------------	-----------------	-------------	----------

nuncTM

Nunc A/S · Kamstrupvej 90 · P.O.Box 280
DK-4000 Roskilde · Denmark
Tel +45 4631 2000 · Fax +45 4631 2175
E-mail: infocity@nunc.dk

www.nuncbrand.com

Apogent

SUBJECT: 110 x $\phi 17$ Disposable Tube.	DRAWN BY/DATE: PS/30.09.96	REVIEW BY/DATE: -	VALID FROM: 30.09.96	APPROVED BY/DATE: JP/30.09.96
	MATERIAL: PS	DRAWING NO.: 2202	VERSION: 2.1	PAGE TYPE: C
SCALE: 2:1				METRIC First Angle Projection.