


SEE BELOW	
WORKING RANGE ( $\mu$ L)	50 – 250 $\mu$ L
MATERIAL	PS
	MM      INCH
A, BASE LENGTH	127,8      5,03
B, BASE WIDTH	85,5      3,35
C, TOP LENGTH	123,3      4,84
D, TOP WIDTH	81,2      3,20
E, FLANGE LONG SIDE	2,4      0,08
F, OVERALL HEIGHT	14,4      0,55
G1, DIA. TOP	$\varnothing$ 7 $\varnothing$ 0,28
G2, DIA. BOTTOM	$\varnothing$ 6,2 $\varnothing$ 0,24
H, A - I LOCATION (X)	14,4      0,55
J, A - I LOCATION (Y)	11,4      0,45
K, WELL TO WELL SPACING	9      0,35
L, WELL DEPTH	11,4      0,45
M, BOTTOM OF WELL DISTANCE	3,2      0,13
N, FLANGE SHORT SIDE	2,4      0,10

CATALOG NO	SURFACE	COLOR	STERILE	WITH LID?	UNITS PK/CAS
436006	AMINO	CLEAR	-	-	5/30
436007	AMINO	WHITE	-	-	5/30
436008	AMINO	BLACK	-	-	5/30
436014	STREPTAVIDIN	CLEAR	-	-	1/15
436015	STREPTAVIDIN	WHITE	-	-	1/15
436016	STREPTAVIDIN	WHITE	-	-	1/15
436024	NICKEL-CHELATE	CLEAR	-	-	1/15
436026	NICKEL-CHELATE	BLACK	-	-	1/15
436027	NICKEL-CHELATE	WHITE	-	-	1/15
436032	GLUTATHIONE	CLEAR	-	-	1/15
436033	GLUTATHIONE	WHITE	-	-	1/15
436034	GLUTATHIONE	BLACK	-	-	1/15

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
  - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



<b>Thermo Fisher</b> <b>SCIENTIFIC</b>	Roskilde Site - Kamstrupvej 90 - P.O.Box 280 DK-4000 Roskilde - Denmark Tlf +45 4631 2000 - Fax +45 4631 2042 E-mail: info.nunc@thermofisher.com			<a href="http://www.nuncbrand.com">www.nuncbrand.com</a>		
	SUBJECT: <b>F96 MICROWELL PLATE</b> <b>IMMOBILIZER, PS</b>	DRAWN BY/DATE:	CHANGE ORDER NO.:	VALID FROM:	APPROVED BY/DATE:	FORMAT:
ITEM NUMBER:	MATERIAL:	DRAWING NO.:	VERSION:	PAGE TYPE:	SCALE:	METRIC 

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC