


CATALOG NUMBER	SEE BELOW	
WORKING RANGE (µL)	10 – 100 µL	
MATERIAL	PS	
	MM	INCH
A. BASE LENGTH	127,6	5,02
B. BASE WIDTH	85,5	3,37
C. TOP LENGTH	123,3	4,85
D. TOP WIDTH	81,2	3,20
E. FLANGE LONG SIDE	2,7	0,11
F. OVERALL HEIGHT	14,4	0,57
G1. WELL WIDTH - TOP	3,7	0,15
G2. WELL WIDTH - BOTTOM	2,9	0,11
H. A - I LOCATION (X)	12,1	0,48
J. A - I LOCATION (Y)	9	0,35
K. WELL TO WELL SPACING	4,5	0,18
L. WELL DEPTH	11,7	0,46
M. BOTTOM OF WELL DISTANCE	2,7	0,11
N. FLANGE SHORT SIDE	2,7	0,11

CATALOG NO.	SURFACE	COLOR	STERILE	WITH LID*	UNITS/PK/CA
436009	AMINO	CLEAR	-	-	5/30
436011	AMINO	WHITE	-	-	5/30
436012	AMINO	BLACK	-	-	5/30
436017	STREPTAVIDIN	CLEAR	-	-	1/15
436018	STREPTAVIDIN	WHITE	-	-	1/15
436019	STREPTAVIDIN	BLACK	-	-	1/15
436028	NICKEL-CHELATE	CLEAR	-	-	1/15
436029	NICKEL-CHELATE	WHITE	-	-	1/15
436031	NICKEL-CHELATE	BLACK	-	-	1/15
436036	GLUTATHIONE	CLEAR	-	-	1/15
436037	GLUTATHIONE	WHITE	-	-	1/15
436038	GLUTATHIONE	BLACK	-	-	1/15

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
  - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



<b>ThermoFisher</b> <b>SCIENTIFIC</b>		Roskilde Site - Kamstrupvej 90 - P.O.Box 280 DK-4000 Roskilde - Denmark Tel +45 4631 2000 - Fax +45 4631 2042 E-mail: Info.nunc@thermofisher.com		<a href="http://www.nuncbrand.com">www.nuncbrand.com</a>	
		SUBJECT: <b>384 WELL PLATE</b> <b>IMMOBILIZER, PS</b>	DRAWN BY/DATE: CHANGE ORDER NO.: VALID FROM:	APPROVED BY/DATE: ITEM NUMBER:	MATERIAL: DRAWING NO.:

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC