



A-A ( 2,5 : 1 )



48,3 ± 0,3



Tolerances unless otherwise specified:  
 Without decimal: ±0.5 mm  
 1 place decimal: ±0.1 mm  
 2 place decimal: ±0.05 mm

**Thermo Fisher**  
 S C I E N T I F I C

Nunc A/S · Kamstrupvej 90 · P.O.Box 280  
 DK-4000 Roskilde · Denmark  
 Tel +45 4631 2000 · Fax +45 4631 2042  
 E-mail: info.nunc@thermofisher.com

[www.thermoscientific.com](http://www.thermoscientific.com)

SUBJECT:

Assembly of 1,0 ml cryo  
 vial & stopper

DRAWN BY/DATE:

-

CHANGE ORDER NO.:

CO-5489

VALID FROM:

07-10-2014

APPROVED BY/DATE:

ROSKGA 07-10-2014

FORMAT:

A4

ITEM NUMBER:

-

MATERIAL:

-

DRAWING NO.:

3558

VERSION:

11-b

PAGE TYPE:

-C

SCALE:

2,5:1

METRIC

First Angle

Projection

