


CATALOG NUMBER	SEE BELOW	
WORKING RANGE (µL)	50 – 250µL	
MATERIAL	PS	
	MM	INCH
A. BASE LENGTH	127,6	5,02
B. BASE WIDTH	85,5	3,37
C. TOP LENGTH	123,5	4,86
D. TOP WIDTH	81,2	3,20
E. FLANGE LONG SIDE	2,7	0,11
F. OVERALL HEIGHT	14,4	0,57
G1, DIA TOP	Ø6,8	Ø0,27
G2, DIA BOTTOM	Ø6,4	Ø0,25
H. A - I LOCATION (X)	14,4	0,57
J. A - I LOCATION (Y)	11,4	0,45
K. WELL TO WELL SPACING	9	0,35
L. WELL DEPTH	9,8	0,39
M. BOTTOM OF WELL DISTANCE	4,6	0,18
N. FLANGE SHORT SIDE	2,7	0,11
*P. LID LENGTH	127,4	5,02
*Q. LID WIDTH	85,4	3,36
*R. LID HEIGHT	9	0,35
S. STACKED HEIGHT – NO LID	27	1,06
*T. STACKED HEIGHT W(LIDS)	31,1	1,22
*U. ASSEMBLY HEIGHT W/LID	16,5	0,65

CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	UNITS/ PK/CASE
277143	NON-TREATED	CLEAR	+	+	1/50
249935	NON-TREATED	CLEAR	+	+	10/80
249940	NON-TREATED	CLEAR	+	+	10/80
249952	NON-TREATED	CLEAR	+	+	10/160
249670	NON-TREATED	CLEAR	-	-	10/180
249652	NON-TREATED	CLEAR	+	-	1/50

\* DIM.S APPLY ONLY WHERE "WITH LID" = +. USING NUNC LID 2205

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
  - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



**ThermoFisher**  
SCIENTIFIC

Roskilde Site • Kamstrupvej 90 • P.O.Box 280  
DK-4000 Roskilde • Denmark  
Tid +45 4631 2000 • Fax +45 4631 2042  
E-mail: info.nunc@thermoFisher.com

[www.nuncbrand.com](http://www.nuncbrand.com)

SUBJECT:	DRAWN BY/DATE:	CHANGE ORDER NO.:	VALID FROM:	APPROVED BY/DATE:	FORM 1:
V96 MICROWELL PLATE, CLEAR, PS					

CONFIDENTIAL  
THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC