




⊙ $\phi 35.8$



⊙ $\phi 34.7$



Tolerances unless otherwise specified:

Without decimal: ± 0.5 mm

1 place decimal: ± 0.1 mm

2 place decimal: ± 0.05 mm

Undefined radius: R 0.20 mm

SUBJECT:

Nundon PetriDish 35 x 10 mm

DRAWN BY/DATE:

MZO/06.12.00

CHANGE ORDER NO.:

CO-0000

VALID FROM:

2007.12.05

APPROVED BY/DATE:

JR2007.11.29

FORMAT:

A3

Roskilde Site • Kamstrupvej 90 • P.O.Box 280
 DK-4000 Roskilde • Denmark
 Tel +45 4631 2000 • Fax +45 4631 2042
 E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

VERSION:

10

DRAWING NO.:

3254

MATERIAL:

PS

PAGE TYPE:

-C

SCALE:

2:1

APPROVED BY/DATE:

JR2007.11.29

FORMAT:

A3

METRIC

ISO

2:1

(2.5:1)