



CATALOG NUMBER	SURFACE	COLOR	STERILE	UNITS PER PACK/CASE
266120	NON-TREATED	CLEAR	YES	25/100

CATALOG NUMBER.....:	SEE ABOVE	
NOM. CAP (µL).....:	400	
MATERIAL.....:	PS	
	INCH	MM
L1 (BASE LENGTH).....:	5.030	127.76
L2 (TOP LENGTH).....:	4.86	123.4
L4 (H-12 LOCATION) (X AXIS).....:	4.46	113.27
L5 (FLANGE LENGTH, LONG SIDE).....:	0.86	21.8
W1 (BASE WIDTH).....:	3.365	85.47
W2 (TOP WIDTH).....:	3.20	81.2
W4 (H-12 LOCATION) (Y AXIS).....:	2.92	74.2
X1 (WELL-TO-WELL SPACING).....:	0.35	9.0
Y1 (WELL-TO-WELL SPACING).....:	0.35	9.0
D1 (TOP OF WELL DIA.).....:	0.28	7.1
D2 (BOTTOM OF WELL DIA.).....:	0.24	6.1
H1 (OVERALL HEIGHT).....:	0.57	14.4
H2 (VERT. DISTANCE FROM INSIDE WELL BTM. TO RESTING PLANE).....:	0.125	3.17
H4 (VERT. DISTANCE FROM OUTSIDE WELL BTM. TO RESTING PLANE).....:	0.085	2.16
H5 (FLANGE HEIGHT, SHORT SIDE).....:	0.09	2.2
H6 (FLANGE HEIGHT, LONG SIDE).....:	0.255	6.5
L3 (A-1 LOCATION) (X AXIS).....:	0.57	14.4
W3 (A-1 LOCATION) (Y AXIS).....:	0.44	11.2
H3 (WELL DEPTH).....:	0.445	11.3
LID L1 (LID LENGTH).....:	5.01	127.35
LID W1 (LID WIDTH).....:	3.36	85.4
LID H1 (EXCLUDING STACK LUGS).....:	0.31	8.0
LID H2 (LID HEIGHT INCLUDES LUGS).....:	0.35	9.0
T UNLIDDED STACK (HEIGHT OF TWO STACKED UNLIDDED PLATES).....:	1.08	27.4
U (HEIGHT OF SINGLE LIDDED PLATE, EXCLUDING STACK LUGS).....:	0.68	17.3
V SINGLE LIDDED PLATE (INCLUDING STACK LUGS).....:	0.64	16.3
X LIDDED STACK (DOES NOT INCLUDE STACK LUG).....:	1.28	32.5
Y LIDDED STACK (INCLUDING LID LUG).....:	1.32	33.5



SPECS ON BACK

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP
CK'D	KOW	11/20/13	

PSS 266120

REV. 1113

DESCRIPTION: 96 WELL, PS FLAT BOTTOM PLATE PINCH BAR