




CATALOG NUMBER	NOMINAL CAPACITY	MATERIAL	
		BOTTLE	CLOSURE
376814	250mL	POLYPROPYLENE	HDPE

NOM BRIM CAPACITY	255 mL	
	INCH	MM
A (OVERALL HEIGHT)	5.68	144.2
B (BOTTLE HEIGHT)	5.48	139.1
C (CLOSURE DIA)	2.38	60.4
D (BOTTLE DIA)	2.36	59.7
E (NECK ID)	1.84	46.8
F (WALL THICKNESS)	.05	1.3

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

ThermoFisher
S C I E N T I F I C

Thermo Fisher Scientific
75 Pondcroft Creek Drive
Rochester, NY 14625-2385, USA

DATE: FMP
TDC: 3/18/15
REV: 0315

PSS 376814

250mL CENTRIFUGE BOTTLE AND CLOSURE ASSEMBLY

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.