




FOR INTERNAL MANUFACTURING PURPOSES ONLY



UNCONTROLLED COPY

BOTTOM RADIUS 8.683mm MAX., REFERENCE

 <b>CHASE SCIENTIFIC GLASS</b> division of CHASE INSTRUMENTS CORP.		
TITLE DISPOSABLE CULTURE TUBE 18 x 125mm	SCALE NOT TO SCALE	DRAWN BY REVISION # 2
DATE 07/19/95	APPROVED BY 	DRAWING NUMBER B1608
GLASS TYPE I - CLEAR - 51 EXP.		