







REV. ECO NO.	DESCRIPTION
A 000024	RELEASE TO PRODUCTION.



COMPATIBLE OF ALL RESERVOIRS.

FOR REFERENCE ONLY		INNOVATIVE CUSTOM MICROPLATES		1998 WESTOVER RD. CHICOPEE, MA 01022 TEL. (413) 593-5300 FAX (413) 593-5000	
© 2001 CONFIDENTIAL TRADE SECRETS PATENTS ISSUED AND PENDING DO NOT COPY OR DISCLOSE WITHOUT WRITTEN PERMISSION OF INNOVATIVE MICROPLATES		TOLERANCES UNLESS OTHERWISE SPECIFIED XX ± .03 XXX ± .015 XXXX ± .005 FRACTION: 1/64 ANGULAR: 30°15'		MATERIAL: RESPPCP NATURAL POLYPROPYLENE CONCENTRATE:	
PRT NAME: RESERVOIR BAFFLE		CATALOG NO.: A50016		DATE: DRAWN BY: APPR BY: PRT SIZE: SCALE: SHEET 1 OF 1 REV. A	

REV. ECO NO.	DESCRIPTION
A 000023	RELEASE TO PRODUCTION.



NOTE:
MAX VOLUME OF RESERVOIR 90 ml

FOR REFERENCE ONLY

INNOVATIVE
CUSTOM MICROPLATES

1998 WESTOVER RD.
CHICOPEE, MA 01022
TEL (413) 593-5300

© 2001
CONFIDENTIAL
TRADE SECRETS
PATENTS ISSUED
AND PENDING
DO NOT COPY
OR DISCLOSE
WITHOUT WRITTEN
PERMISSION

TOLERANCES UNLESS
OTHERWISE SPECIFIED
XX ± .03
XXX ± .015
XXXX ± .0005
FRACTION: 1/64
ANGULAR: 10° 15'
RADIUS 0.015" IN LESS

MATERIAL: RESPPCP
NATURAL POLYPROPYLENE
CONCENTRATE:

PRT NAME:
RESERVOIR, LOW PROFILE, FLAT BOTTOM, NATURAL PP

CATALOG NO.: S30023