



REV.	ECO NO.	DESCRIPTION
A	000092	RELEASE TO PRODUCTION.



NOTE:  
1) MAX VOLUME OF COLUMN 7ml

 <b>Seahorse Bioscience</b> 300 Griffith Road Chicopee, MA 01022 Tel: (413) 593-5300 Fax: (413) 593-5900	© 2001 CONFIDENTIAL TRADE SECRETS PATENTS ISSUED AND PENDING DO NOT COPY OR DISCLOSE WITHOUT WRITTEN PERMISSION OF INNOVATIVE MICROPLATE	TOLERANCES UNLESS OTHERWISE SPECIFIED XX ± .03 XXX ± .015 XXXX ± .0005 FRACTION: ± 1/64 ANGULAR: 40°/15' RADIUS 0.015" UNLESS OTHERWISE SPECIFIED BREAK ALL SHARP EDGES	<b>MATERIAL:</b> RESPPCP NATURAL POLYPROPYLENE <b>CONCENTRATE:</b>	<b>PRT NAME:</b> 12 COLUMN, LOW PROFILE RESERVOIR, NATURAL PP	<b>PRT NO.:</b> S30028	<b>SCALE:</b> NTS	<b>SHEET</b> 1 OF 1	<b>REV.</b> A
	<b>DATE:</b> 1/3/03 <b>DRAWN BY:</b> S. MORLEY <b>APPR BY:</b> S. YOUNG	<b>PRT SIZE:</b> A	<b>SCALE:</b> NTS	<b>SHEET</b> 1 OF 1	<b>REV.</b> A			