


REV.	ECO NO.	DESCRIPTION
A	000316	RELEASE TO PRODUCTION.
B	000325	ADDED CENTER WALL



NOTE:  
 1.) VOLUME RESERVOIR 155ml  
 2.) VOLUME CONTROLLED WELL 2.60ml

FOR REFERENCE ONLY

**Seahorse Bioscience**  
 300 Griffith Road  
 Chicopee, MA 01022  
 Tel: (413) 593-5300  
 Fax: (413) 593-5900

© 2001 CONFIDENTIAL TRADE SECRETS PATENTS ISSUED AND PENDING DO NOT COPY OR DISCLOSE WITHOUT WRITTEN PERMISSION OF INNOVATIVE MICROPLATE	TOLERANCES UNLESS OTHERWISE SPECIFIED XX ± .03 XXX ± .015 XXXX ± .0005 FRACTION: ± 1/64 ANGULAR: ±0°15' RADIUS 0.015" UNLESS OTHERWISE SPECIFIED BREAK ALL SHARP EDGES	<b>MATERIAL:</b> RESPPCP NATURAL POLYPROPYLENE		<b>PRT NAME:</b> 384 WELL RESERVOIR WITH FOUR CONTROLLED COLUMNS	
		<b>CONCENTRATE:</b>		<b>CATALOG NO.:</b> S30041	
<b>DATE:</b> 3/3/05		<b>DRAWN BY:</b> KN		<b>APPR BY:</b>	
<b>PRT SIZE:</b> A		<b>SCALE:</b> NTS		<b>SHEET</b> 1 of 1	
<b>REV.</b> B					