


Component Schematic

384-well Microplate (PP) Deep Well



L	128.15 ± 0.10	Hr	6.75 ± 0.10
W	85.70 ± 0.10	Wr	21.30 ± 0.10
H	30.20 ± 0.05	P	4.50 ± 0.05
Hs	2.55 ± 0.10	D	3.65 ± 0.05
X1	12.10 ± 0.05	Hw	26.35 ± 0.05
Xn	115.60 ± 0.05	C	24
Y1	9.10 ± 0.05	R	16
Yn	76.60 ± 0.05		

Curvature, relative to resting plane ≤0.20

Volume	300µl
Well form	V Bottom
Material	PP