

ANIOS'CLEAN EXCEL D

Stabilité et conditions de conservation

Stability and storage conditions

Produit pur :

- Conservation entre +5°C et +25°C
- Stabilité :

- Bidon et flacon doseur non ouvert : 24 mois à partir de la date de fabrication indiquée sur l'étiquette.
- Bidon et flacon doseur ouvert : maintien de la durée de conservation de 24 mois à condition de refermer le flacon après utilisation.
- Doses : 18 mois à partir de la date de fabrication indiquée sur l'étiquette.

Produit dilué:

Remplacer la solution dès que celle-ci est visiblement souillée et auminimum une fois par jour pour maintenir l'efficacité du produit.

Pure product:

- Storage between +5°C and +25°C
- Stability :

- Not opened dosing bottle and can : 24 months from the manufacturing date indicated on the label.
- Opened dosing bottle and can : maintenance of the 24 months in condition if closing the bottle after use.
- Dose: 18 months from the conditioning date indicated on the label.

Diluted product :

Replace the solution as soon as it is visibly soilde and at least once per day to maintain the effectiveness of the product.